

Human KLRC1 & CD94 Heterodimer protein (Recombinant) (N-8*His & N-Flag) (STJP013271)

GENERAL INFORMATION

| | |
|---------------|-----------------|
| Product Type | Proteins |
| Applications | ELISA/WB |
| Host / Source | Mammalian Cells |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4. |
| Molecular Weight | 34.4 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|------------------------------------|
| Immunogen Region | [Arg100-Leu233] & [Ser34-Ile179] |
|------------------|------------------------------------|

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|