

## Human CD123/IL3RA protein (Recombinant) (C-Fc) (STJP013275)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilized from a 0.22 µm filtered solution in PBS, pH 7.4.  |
| Molecular Weight    | 60 kDa  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">3563</a>        |
| Gene Symbol      | <a href="#">IL3RA</a>       |
| UniProt ID       | <a href="#">IL3RA_HUMAN</a> |
| Immunogen Region | Thr19-Arg305                |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|