

## Human CST5 protein (Recombinant) (C-His) (STJP013281)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 6% Trehalose, 2% Mannitol, 0.05% Tween 80, pH 8.5.                                                 |
| Molecular Weight    | 14.9 kDa                                                                                                                                                            |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">1473</a>       |
| Gene Symbol      | <a href="#">CST5</a>       |
| UniProt ID       | <a href="#">CYTD_HUMAN</a> |
| Immunogen Region | Gly21-Val142               |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|