

Human EPHA8 protein (Recombinant) (C-Fc) (STJP013290)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, pH 7.4.
Molecular Weight	83.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2046
Gene Symbol	EPHA8
UniProt ID	EPHA8_HUMAN
Immunogen Region	Glu31-Thr542

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--