

Human CD362/SDC2 protein (Recombinant) (C-His) (STJP013309)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	16.86 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 °C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	6383
Gene Symbol	SDC2
UniProt ID	SDC2_HUMAN
Immunogen Region	Met1-Glu144

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--