

Human CD147/BSG protein (Recombinant) (C-mFc) (STJP013317)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 20mM PB, 150mM NaCl, pH 7.2.
Molecular Weight	46.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Ala22-His205
------------------	--------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--