

Human WNT3A protein (Recombinant) (N-fusion design) (STJP013383)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 10mM PB, 5% Sucrose, 4% Mannitol, 0.01% Tween 80, pH 7.4.
Molecular Weight	105.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	89780
Gene Symbol	WNT3A
UniProt ID	WNT3A_HUMAN
Immunogen Region	Ser19-Lys352

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--