

## Human BTN1A1 protein (Recombinant) (C-Fc) (STJP013444)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4.                                                                                                      |
| Molecular Weight    | 50.8 kDa                                                                                                                                                            |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">696</a>         |
| Gene Symbol      | <a href="#">BTN1A1</a>      |
| UniProt ID       | <a href="#">BT1A1_HUMAN</a> |
| Immunogen Region | Ala27-Arg242                |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|