

Human ACVR2B/ACTR-IIB protein (Recombinant) (C-Fc-His) (STJP013456)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 10% Trehalose, 3% Mannitol, 0.05% Tween 80, 10 mM Methionine, pH 8.5.
Molecular Weight	41.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	93
Gene Symbol	ACVR2B
UniProt ID	AVR2B_HUMAN
Immunogen Region	Ser19-Thr134

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--