

Human PDIA6 protein (Recombinant) (C-His) (STJP013476)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Molecular Weight	47.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	10130
Gene Symbol	PDIA6
UniProt ID	PDIA6_HUMAN
Immunogen Region	Leu20-Leu440
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.