

Human CD279/PDCD1/PD1 protein (Recombinant) (C-hFc-His) (STJP013478)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 5% Trehalose, 5% Mannitol.
Molecular Weight	42.80 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	5133
Gene Symbol	PDCD1
UniProt ID	PDCD1_HUMAN
Immunogen Region	Leu25-Thr168

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--