

Human ANGPTL8 protein (Recombinant) (N-Fc) (STJP013510)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 10% Trehalose, 0.05% Tween 80, pH 8.0.
Molecular Weight	46 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	55908
Gene Symbol	ANGPTL8
UniProt ID	ANGL8_HUMAN
Immunogen Region	Ala22-Ala198

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--