

Human CCL17/TARC protein (Recombinant) (C-His) (STJP013541)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in 20mM His-HCl, 10% Trehalose, 5% Maltose, 100mM NaCl, 0.1% Tween 80, pH 6.5.
Molecular Weight	9.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	6361
Gene Symbol	CCL17
UniProt ID	CCL17_HUMAN
Immunogen Region	Ala24Ser94

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------