

## Human CD353/SLAMF8 protein (Recombinant) (C-His) (STJP013601)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	26.87 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 °C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">56833</a>
Gene Symbol	<a href="#">SLAMF8</a>
UniProt ID	<a href="#">SLAF8_HUMAN</a>
Immunogen Region	Met1-Asp233

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--