

## Human CD57/B3GAT1 protein (Recombinant) (N-His) (STJP013603)

### GENERAL INFORMATION

|               |                 |
|---------------|-----------------|
| Product Type  | Proteins        |
| Applications  | ELISA/WB        |
| Host / Source | Mammalian Cells |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilized from a 0.22 µm filtered solution in 20mM Citrate, 8% Sucrose, 100mM NaCl, 0.05% Tween 80, pH 6.0.   |
| Molecular Weight    | 36.2 kDa  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">27087</a>       |
| Gene Symbol      | <a href="#">B3GAT1</a>      |
| UniProt ID       | <a href="#">B3GA1_HUMAN</a> |
| Immunogen Region | His25-Ile334                |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|