

Human DLL3 protein (Recombinant) (C-His) (STJP013639)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 100mM NaCl, 8% Sucrose, 20% Glycerol, 0.05% Tween 80, pH 6.0.
Molecular Weight	15.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10683
Gene Symbol	DLL3
UniProt ID	DLL3_HUMAN
Immunogen Region	Arg353-Leu492

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--