

Human SIGLEC15/CD33L3 protein (Recombinant) (C-mFc) (STJP013642)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in PBS, 150mM NaCl, 0.3% Chaps, 5% Trehalose, pH 7.4.
Molecular Weight	52.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	284266
Gene Symbol	SIGLEC15
UniProt ID	SIG15_HUMAN
Immunogen Region	Phe20-Thr263

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--