

Human HLA-A*0201 HBs Ag Complex protein (Recombinant-Active) (C-10*His) (STJP013676)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 500mM NaCl, 5% Trehalose, 50% Glycerol, 0.06% Tween 80, pH 7.4.
Molecular Weight	49 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	[ILSPFLPLL] & [Ile21-Met119] & [Gly25-Ile308 (A269V)]
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--