

Human GALNT7 protein (Recombinant) (C-His) (STJP013720)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM Tris-HCl, 10mM reduced Glutathione, pH 8.0.
Molecular Weight	73.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	51809
Gene Symbol	GALNT7
UniProt ID	GALT7_HUMAN
Immunogen Region	Pro30-Val657

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--