

Mouse LDHA protein (Recombinant) (C-His) (STJP013780)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 50% Glycerol, 500 mMNaCl, 5% Trehalose, 5% Mannitol, 0.01% Tween 80, 1mM EDTA, pH9.0.
Molecular Weight	37.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	16828
Gene Symbol	Ldha
UniProt ID	LDHA_MOUSE
Immunogen Region	Met1-Phe332

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--