

Mouse Latent TGF-beta-1 protein (Recombinant) (N-8*His) (STJP013786)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 15% Trehalose, 4% Mannitol, 0.1% Tween 80, pH 8.2.
Molecular Weight	12.8 kDa & 31.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	21803
Gene Symbol	Tgfb1
UniProt ID	TGFB1_MOUSE
Immunogen Region	Leu30-Ser390 (C33S)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--