

Mouse CD142/F3/TF protein (Recombinant) (C-His) (STJP013838)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 150mM NaCl, 0.05% Brij35, pH 7.5.
Molecular Weight	26.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	14066
Gene Symbol	F3
UniProt ID	TF_MOUSE
Immunogen Region	Ala29-Glu251
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.