

## Mouse CD79 Beta/CD79B protein (Recombinant) (C-Fc) (STJP013915)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in PBS, pH 7.4.
Molecular Weight	42 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">15985</a>
Gene Symbol	<a href="#">Cd79b</a>
UniProt ID	<a href="#">CD79B_MOUSE</a>
Immunogen Region	Val26-Asp158

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--