

Mouse SCARB2/CD36L2 protein (Recombinant) (C-Fc) (STJP013961)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 50mM Tris-Citrate, 0.3M NaCl, pH 6.5.
Molecular Weight	73.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	12492
Gene Symbol	Scarb2
UniProt ID	SCRB2_MOUSE
Immunogen Region	Arg27-Thr432

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--