

Mouse ACE2 protein (Recombinant) (C-10*His) (STJP014006)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 25mM Tris-HCl, 150mM NaCl, 1mM ZnCl ₂ , pH 7.5.
Molecular Weight	84.9 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	70008
Gene Symbol	Ace2
UniProt ID	ACE2_MOUSE
Immunogen Region	Gln18-Thr740

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--