

Mouse ENTPD2/CD39L1 protein (Recombinant) (C-His) (STJP014045)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian Cells
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM Tris-HCl, 10mM CaCl ₂ , 150mM NaCl, 10% Glycerol, pH 7.5.
Molecular Weight	49.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	12496
Gene Symbol	Entpd2
UniProt ID	ENTP2_MOUSE
Immunogen Region	Cys26-Ser462
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.