

Cynomolgus monkey CD40/TNFRSF5 protein (Recombinant) (C-Fc) (STJP014131)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/WB
	Host / Source	Mammalian Cells
PRODUCT PROPERTIES		
	Formulation	Lyophilized from a 0.22 Mu m filtered solution in 50mM Tris-HCl, 100mM Glycine, pH 7.5.
	Molecular Weight	46.3 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Immunogen Region	Glu21-Arg193
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.