

Cynomolgus monkey GDF15 protein (Recombinant) (N-8*His-Flag) (STJP014164)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	Mammalian cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 Mu m filtered solution of 20mM Glycine-HCl, 10% Trehalose, 0.02% Tween 80, pH 3.0.
Molecular Weight	14.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Ala197-Val308
------------------	---------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--