

## Biotin conjugated Cyclic Citrullinated Peptide/CCP protein {Biotin} (STJP014242)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA

### PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in PBS pH 7.4, 4% Trehalose, 1% Mannitol.
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--