

Human PDCD5 protein (Recombinant) (N-His) (STJP014255)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Tris-HCl, 5% Sucrose, 5% Mannitol, 50mM NaCl, 0.05% Tween 80, pH 8.4.
Molecular Weight	16.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	9141
Gene Symbol	PDCD5
UniProt ID	PDCD5_HUMAN
Immunogen Region	Met1-Tyr125

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--