

Human BRD4 protein (Recombinant) (N-10*His-Flag) (STJP014286)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM HEPES, 200mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.
Molecular Weight	49 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	23476
Gene Symbol	BRD4
UniProt ID	BRD4_HUMAN
Immunogen Region	Glu49-Glu460
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.