

Human IDH1 (R132H) protein (Recombinant) (C-8*His) (STJP014294)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 50mM Tris-HCl, 4% Sucrose, 50% Glycerol, 0.02% Tween 80, pH 8.0.
Molecular Weight	48.1 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3417
Gene Symbol	IDH1
UniProt ID	IDHC_HUMAN
Immunogen Region	Met1-Leu414 (R132H)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--