

Human BCAS2 protein (Recombinant) (N-T7 & C-His) (STJP014298)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 20mM Tris-HCl, 200mM NaCl, 2mM DTT, pH 8.0.
Molecular Weight	28.58 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10286
Gene Symbol	BCAS2
UniProt ID	SPF27_HUMAN
Immunogen Region	Ala2-Phe225

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--