

Human AIFM1/AIF protein (Recombinant) (N-His) (STJP014311)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, 50% Glycerol, 2mM EDTA, 0.5M Argine, 5% Trehalose, pH 7.4.
Molecular Weight	56.2 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	9131
Gene Symbol	AIFM1
UniProt ID	AIFM1_HUMAN
Immunogen Region	Glu121-Asp613

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--