

## Human CLIC3 protein (Recombinant) (C-His) (STJP014312)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 200mM Sucrose, 0.05% Tween 80 (w/v) , pH 8.5.                                                            |
| Molecular Weight    | 27.7 kDa                                                                                                                                                              |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">9022</a>        |
| Gene Symbol      | <a href="#">CLIC3</a>       |
| UniProt ID       | <a href="#">CLIC3_HUMAN</a> |
| Immunogen Region | Met1-Arg236                 |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|