

Human CA1/Carbonic anhydrase 1 protein (Recombinant) (C-His) (STJP014318)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 12.5mM Tris-HCl, 75mM NaCl, pH 7.5.
Molecular Weight	29.93 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	759
Gene Symbol	CA1
UniProt ID	CAH1_HUMAN
Immunogen Region	Ala2-Phe261

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--