

Human ASS1 protein (Recombinant) (N-His) (STJP014319)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, 50mM Imidazole, 1mM DTT, 40% Glycerol, pH 7.5.
Molecular Weight	42.8 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	445
Gene Symbol	ASS1
UniProt ID	ASSY_HUMAN
Immunogen Region	Met1-Lys412

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--