

## Human CALCA/Calcitonin protein (Recombinant) (C-His) (STJP014340)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Supplied as a 0.22 µg/ml filtered solution in 20mM PB, 150mM NaCl, 50% Glycerol, pH 7.4.  |
| Molecular Weight    | 10.5 kDa  |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">796</a>        |
| Gene Symbol      | <a href="#">CALCA</a>      |
| UniProt ID       | <a href="#">CALC_HUMAN</a> |
| Immunogen Region | Tyr58-Asn141               |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|