

Human DNTT protein (Recombinant) (N-His) (STJP014359)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 0.1M KH ₂ PO ₄ , 200mM KCl, 1mM 2-Mercaptoethanol, 50% Glycerol, pH 7.2.
Molecular Weight	60.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	1791
Gene Symbol	DNTT
UniProt ID	TDT_HUMAN
Immunogen Region	Met1-Ala509

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--