

Human TUBB4A protein (Recombinant) (N-His) (STJP014364)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Histidine-HCl, 10% Trehalose, 100mM NaCl, 0.05% Tween 80, pH 6.5.
Molecular Weight	51.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	10382
Gene Symbol	TUBB4A
UniProt ID	TBB4A_HUMAN
Immunogen Region	Met1-Ala444

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--