

## Human HBQ1 protein (Recombinant) (N-His) (STJP014385)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 20mM Histidine, 10% Sucrose, 3% Mannitol, 0.5mM EDTA, 0.05% Tween 80, pH 5.5.
Molecular Weight	17.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3049</a>
Gene Symbol	<a href="#">HBQ1</a>
UniProt ID	<a href="#">HBAT_HUMAN</a>
Immunogen Region	Met1-Arg142

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--