

Human Gal1/LGALS1 protein (Recombinant) (C-His) (STJP014387)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mu m filtered solution in 20mM Citrate, 8% Trehalose, 4% Mannitol, 50mM NaCl, 1mM TCEP, 0.05% Tween 80, pH 4.0.
Molecular Weight	15.78 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3956
Gene Symbol	LGALS1
UniProt ID	LEG1_HUMAN
Immunogen Region	Ala2-Asp135

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--