

Human HMGB1 protein (Recombinant) (Tag free (N-terminal)) (STJP014388)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 Mg/ml filtered solution in 50mM HEPES, 500mM NaCl, 0.5mM DTT, pH 7.9.
Molecular Weight	10.4 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	3146
Gene Symbol	HMGB1
UniProt ID	HMGB1_HUMAN
Immunogen Region	Met1-Phe89

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--