

## Human HMGB1 protein (Recombinant) (Tag free) (STJP014389)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 M filtered solution in 50mM HEPES-Na, 500mM NaCl, 0.6mM DTT, pH 7.9.
Molecular Weight	10 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3146</a>
Gene Symbol	<a href="#">HMGB1</a>
UniProt ID	<a href="#">HMGB1_HUMAN</a>
Immunogen Region	Pro92-Val176

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--