

Human SULT1A3 protein (Recombinant) (N-His) (STJP014396)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 6% Sucrose, 4% Mannitol, 0.05% Tween 80, pH 8.5.
Molecular Weight	35.6 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	445329/6818
Gene Symbol	SULT1A3
UniProt ID	ST1A3_HUMAN
Immunogen Region	Met1-Leu295

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--