

Human ODC1 protein (Recombinant) (N-T7 & C-His) (STJP014408)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris-HCl, 100mM NaCl, 8% Trehalose, 2% Mannitol, 4 mM TCEP, 0.05% Tween 80, pH 8.0.
Molecular Weight	53.5 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	4953
Gene Symbol	ODC1
UniProt ID	DCOR_HUMAN
Immunogen Region	Met1-Val461
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.