

Human FABP2 protein (Recombinant) (N-His & C-His) (STJP014409)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in 20mM PB, 50mM NaCl, 8% Trehalose, 0.05% Tween 80, pH 6.0.
Molecular Weight	18.44 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	2169
Gene Symbol	FABP2
UniProt ID	FABPI_HUMAN
Immunogen Region	Met1-Asp132

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--