

Human CCL2/MCP-1 protein (Recombinant) (Tag free) (STJP014413)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from a 0.22 µm filtered solution in 20mM PB, 6% Sucrose, 4% Mannitol, 0.05% Tween 80, pH 6.0.
Molecular Weight	8.7 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	6347
Gene Symbol	CCL2
UniProt ID	CCL2_HUMAN
Immunogen Region	Gln24-Thr99

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--