

## Human FABP4 protein (Recombinant) (N-His) (STJP014421)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | ELISA/WB |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilized from a 0.22 Mu m filtered solution in PBS, pH 7.4.                                                                                                      |
| Molecular Weight    | 16.88 kDa                                                                                                                                                           |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at-20 to-80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">2167</a>        |
| Gene Symbol      | <a href="#">FABP4</a>       |
| UniProt ID       | <a href="#">FABP4_HUMAN</a> |
| Immunogen Region | Cys2-Ala132                 |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|