

Human RO52 protein (Recombinant) (N-His) (STJP014438)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as a 0.2 µg/ml filtered solution of 20mM Tris-HCl, 10% Trehalose, 100mM Arg-HCl, 0.05% Tween 80, pH8.5.
Molecular Weight	56.3 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	6737
Gene Symbol	TRIM21
UniProt ID	R052_HUMAN
Immunogen Region	Met1-Tyr475

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--